









Ameda Replacement Parts

Ameda is considered by many Canadian health care professionals to be the premier breast pump manufacturer in the world today. With over 70 years of engineering know-how and proven success, millions of new mothers and their babies have benefited from Ameda breastfeeding pumps.



Replacement Parts and Accessories

Description	
	<p>CustomFit™ Breast Flange M-L (28.5mm / 30.5mm) Contains 2 larger breast flanges and special reducing inserts to create a smaller size for maximum pumping effectiveness. The 28.5 mm insert fits inside the 30.5 mm flange. No special connectors needed.</p> <p><i>Includes: (2) 30.5mm Breast Flanges and (2) 28.5mm Inserts.</i></p>
	<p>CustomFit™ Breast Flange XI-XXL (32.5mm / 36mm) Contains 2 larger breast flanges and special reducing inserts to create a smaller size for maximum pumping effectiveness. The 32.5 mm insert fits inside the 36.0 mm flange. No special connectors needed.</p> <p><i>Includes: (2) 36mm Breast Flanges and (2) 32.5mm Inserts.</i></p>
	<p>Flexishield™ Areola Stimulator (2 Pack) Soft, silicone insert fits the standard size Ameda breast flange. Reduces flange size to 21.0 mm for maximum effectiveness with smaller nipples. Helps stimulate the let-down reflex. Works with any 25 mm breast flange.</p> <p><i>Includes: (2) Flexishields.</i></p>
	<p>Reducing Insert (22.5mm) Reduces flange size to 22.5 mm for maximum effectiveness for smaller nipples. Fits any 25 mm flange. (not in retail packaging).</p> <p><i>Includes: (1) Reducing Insert.</i></p>
	<p>Spare Parts Kit (fits all Ameda HygieniKits)</p> <p><i>Includes: (4) Valves, (2) Silicone Tubing, (2) Silicone Diaphragms, (2) Adapter Caps, (1) Tubing Adapter.</i></p>
	<p>Valves (2 Pack) 2 HygieniKit replacement valves. Each HygieniKit comes with a spare set. Valves should be replaced every 6 weeks depending on use.</p> <p><i>Includes: (2) Replacement valves.</i></p>

Replacement Parts and Accessories

Description	
	<p>Diaphragms (2 Pack) 2 HygieniKit replacement diaphragms.</p>
	<p>Dual Connector 1 Dual Connector (not in retail packaging).</p>
	<p>Replacement Flange (25mm) Replacement for standard size flange (not in retail packaging).</p>
	<p>Adaptor Cap 1 adaptor cap (not in retail packaging).</p>
	<p>Tubing 1 piece of tubing.</p>
	<p>Car Adapter - Purely Yours 12V adapter can be used to power all Purely Yours and Purely Yours Ultra models.</p>
	<p>Car Adapter - Finesse 9V adapter can be used to power all Finesse models.</p>
	<p>AC Adapter Designed for use with the Ameda Finesse, Purely Yours and Purely Yours Ultra Electric Breast Pumps. <i>Other models available for use in the following areas: North America, United Kingdom, Europe & Australia.</i></p>