














Replacement Parts and Accessories

Ameda is considered by many Canadian health care professionals to be the premier breast pump manufacturer in the world today. With over 70 years of engineering know-how and proven success, millions of new mothers and their babies have benefited from Ameda breastfeeding pumps.

DESCRIPTION	
	<p>CustomFit™ Breast Flange M-L (28.5mm / 30.5mm) Contains 2 larger breast flanges and special reducing inserts to create a smaller size for maximum pumping effectiveness. The 28.5 mm insert fits inside the 30.5 mm flange. No special connectors needed.</p> <p><i>Includes: (2) 30.5mm Breast Flanges and (2) 28.5mm Inserts.</i></p>
	<p>CustomFit™ Breast Flange XI-XXL (32.5mm / 36mm) Contains 2 larger breast flanges and special reducing inserts to create a smaller size for maximum pumping effectiveness. The 32.5 mm insert fits inside the 36.0 mm flange. No special connectors needed.</p> <p><i>Includes: (2) 36mm Breast Flanges and (2) 32.5mm Inserts.</i></p>
	<p>Flexishield™ Areola Stimulator (2 Pack) Soft, silicone insert fits the standard size Ameda breast flange. Reduces flange size to 21.0 mm for maximum effectiveness with smaller nipples. Helps stimulate the let-down reflex. Works with any 25 mm breast flange.</p> <p><i>Includes: (2) Flexishields.</i></p>
	<p>Reducing Insert (22.5mm) Reduces flange size to 22.5 mm for maximum effectiveness for smaller nipples. Fits any 25 mm flange. (not in retail packaging).</p> <p><i>Includes: (1) Reducing Insert.</i></p>
	<p>Spare Parts Kit (fits all Ameda HygieniKits)</p> <p><i>Includes: (4) Valves, (2) Silicone Tubing, (2) Silicone Diaphragms, (2) Adapter Caps, (1) Tubing Adapter.</i></p>
	<p>Valves (2 Pack) 2 HygieniKit replacement valves. Each HygieniKit comes with a spare set. Valves should be replaced every 6 weeks depending on use.</p> <p><i>Includes: (2) Replacement valves.</i></p>

Replacement Parts and Accessories

DESCRIPTION	
	<p>Diaphragms (2 Pack) 2 HygieniKit replacement diaphragms.</p>
	<p>Dual Connector 1 Dual Connector (not in retail packaging).</p>
	<p>Replacement Flange (25mm) Replacement for standard size flange (not in retail packaging).</p>
	<p>Adaptor Cap 1 adaptor cap (not in retail packaging).</p>
	<p>Tubing 1 piece of tubing.</p>
	<p>Car Adapter - Purely Yours 12V adapter can be used to power all Purely Yours and Purely Yours Ultra models.</p>
	<p>Car Adapter - Finesse 9V adapter can be used to power all Finesse models.</p>
	<p>AC Adapter Designed for use with the Ameda Finesse, Purely Yours and Purely Yours Ultra Electric Breast Pumps. <i>Other models available for use in the following areas: North America, United Kingdom, Europe & Australia.</i></p>